

TECHLEM CADAVER STRETCHER



www.midmed.com.au



1300 643 633

TECHLEM CADAVER STRETCHER

Featuring flap access for personal belongings, and an easy-clean base with integrated storage, our cadaver stretcher is a simple, functional design for use in cadaver handling up to 600lbs (c) or 750lbs (brc)

Features:

- Mint Green Vinyl Covered Canopy With
- Flap Access for Personal Belongings
- Light Weight Removable Canopy Frame
- Stainless Steel Recessed Litter
- Techlem “Posi-Lock” Braking System
- 8” (20cm) All Surface, 4 Wheel Central Lock Casters
- “Accu-Steer” Fifth Wheel Steering
- 4 Corner Access Brake / Steer Pedals
- Raise / Lower Pedals – 4 Side Access
- Mint Green Wraparound Bumper & Accents
- Easy Clean Base with Integrated Storage



DIMENSIONS (5000C)

TOTAL DIMENSIONS: 78" X 29"
RECESSED LITTER: 77" X 29"
CANOPY: 76" X 28" X 14"
WEIGHT CAPACITY: 700LBS / 318KG

DIMENSIONS (5000BRC)

TOTAL DIMENSIONS: 81" X 32"
RECESSED LITTER: 80" X 32"
CANOPY: 80" X 32" X 24"
WEIGHT CAPACITY: 750LBS / 341KG

MOVEMENT

HEIGHT: 22" TO 35^{1/2}"(C) OR 22" TO 32" (BRC)
ACCU-STEER 5TH WHEEL
POSI-LOCK BRAKING SYSTEM
8" ALL-SURFACE, LOCKING CASTERS
OPTIONAL FOUR-SIDED ACCESS

SPECIAL FEATURES

-VINYL CANOPY COVER WITH FLAP ACCESS
-LIGHTWEIGHT, REMOVABLE CANOPY FRAME
-STAINLESS STEEL RECESSED LITTER SHEET