

WRASFER™

The new 2 in 1 Disposable Patient Wrap & Transfer sheet

A hygienic and safe solution in patient bedding and transfer

Easy & Safe patient transfer
Cross infection control
Time & Cost saving



WRASFER™


The new 2 in 1 Disposable Patient Wrap & Transfer sheet

A hygienic and safe solution in patient bedding and transfer

R.T.S. disposable patient wrap & transfer sheet is designed for use at paramedic services, hospital emergency wards and nursing homes, in situations where clean emergency bedding is required and patients need to be transported and transferred. Through its integrated fluff pad fluids up to a volume of 5L can be reliably absorbed and contained and patients up to a max. weight of 200 KG can be easily and safely transferred from one carrier to another. It addresses two critical issues: contamination control and easy patient transfer.

The benefits of WRASFER Tm in overview

- Super strong carrying strength up to 200 KG (tested and approved)
- High absorbency of fluids up to 5L (tested and approved)
- Same carrying strength after absorbing 5L fluids (tested and approved)
- Creating a hygienic ambulance and working environment
- Minimizing the risk of cross infection for patients and staff
- Significantly reduced cleaning and turnaround time
- Individually wrapped packs with instruction for use
- Highly versatile application fields
- Very competitive product costs



Product data sheet



Main sheet 220 x 110 CM
fluff pad 200 x 50 CM



3 layers with PE-Foil, nonwoven and cellulose fluff



350 g



Lifting weight up to 200KG



Absorbency up to 5L






Blue-White



Non sterile



Individually packed, 20 in a carton box



exclusive distributors



Unit 4, 62 Borthwick Avenue,
Murarrie, 4172, QLD

PO Box 508
MORNINGSIDE, QLD 4170