

ENFit™ Enteral Feeding Catalogue

2026 EDITION




SALES@MIDMED.COM.AU

The History of ENFit™



In 2006, following several incidences of tubing misconnections within the clinical setting the International Organisation, GEDSA was formed. After years of research and development, a unique locking system was developed called ENFit™. The recommendations for the system were ISO:80369 and this has now evolved to the standard ISO:80369-3.

GEDSA is a non-profit trade association, formed to establish a voice for addressing issues that faced enteral device manufacturers, suppliers and distributors, and to help international standards for healthcare tubing connectors comprise of manufacturers, clinical and medical staff and suppliers worldwide, GEDSA facilitates information flow about the initiative which is designed to increase patient safety and optimal delivery of enteral feeding by reducing the risk of tubing misconnections. Our staff are GEDSA trained to provide you with information in accordance with ISO:80369-3



ENFit™ SINGLE AND HOME USE SYRINGES



STERILE ENFit™ HOME USE SYRINGES - Up to 14 Day Use

Code	Size	Units/Box
MS.SE01H	1ml	100
MS.SE03H	3ml	100
MS.SE05H	5ml	100
MS.SE10H	10ml	100
MS.SE20H	20ml	50
MS.SE60H	60ml	30



- ISO 80369-3 standard tip for enteral use connection
- Silicone O-ring gasket design for clean and multiple use
- Available for Low Dose Tip design on 1ml, 3ml, 5ml
- Purple plunger for ENFit™ enteral identification
- Permanent and heat-etched dual graduations increase safety
- Safe plunger backstop to prevent withdrawal
- Sterilised by EO or Gamma, non-toxic, non-pyrogenic

STERILE ENFit™ SINGLE-USE SYRINGES

Code	Size	Units/Box
MS.SE01	1ml	100
MS.SE03	3ml	100
MS.SE05	5ml	100
MS.SE10	10ml	100
MS.SE20	20ml	50
MS.SE60	60ml	30



- ISO 80369-3 standard tip for enteral use connection
- Available for Low Dose Tip design on 1ml, 3ml and 5ml
- Purple plunger for ENFit™ enteral identification
- Permanent and heat-etched dual graduations increase safety
- Safe plunger backstop to prevent withdrawal
- Sterilised by EO or Gamma, non-toxic, non-pyrogenic

SYRINGE CAP PUSH FIT ENFit™/ORAL



The ENFit™ Syringe Cap Push Fit is designed to provide a secure and reliable closure for both ENFit™ and oral syringes. This versatile cap ensures the safe storage and transportation of syringes, preventing contamination and maintaining the integrity of the contents.

- **Compatibility:** Suitable for both ENFit™ and oral syringes, offering flexibility in usage.
- **Secure Fit:** Push-fit design ensures a tight seal, preventing leaks and spills.
- **Safety:** Made from high-quality, medical-grade materials that are DEHP, latex, and BPA-free.
- **Ease of Use:** Simple to attach and remove, making it convenient for healthcare professionals
- **Hygienic:** Helps maintain sterility of syringe contents, reducing the risk of contamination.

Stock Code: ENT.PFC-100



SHORT TERM PVC FEEDING TUBE (7-10 Days)




Product Features

- Range of Fr sizes (4Fr - 14 Fr)
- Range of lengths (50cm - 92cm)
- 2 non-opposing outlet holes
- Fully radiopaque
- Made of PVC
- Use for up to 10 days and 7 days for neonates
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA Free
- 25 units per box

Code	Fr/CH	Length
ENT.CF01.04.050ISOSAF	4Fr	50cm
ENT.CF01.04.075ISOSAF	4Fr	75cm
ENT.CF01.05.050ISOSAF	5Fr	50cm
ENT.CF01.05.075ISOSAF	5Fr	75cm
ENT.CF01.06.050ISOSAF	6Fr	50cm
ENT.CF01.06.075ISOSAF	6Fr	75cm
ENT.CF01.08.050ISOSAF	8Fr	50cm
ENT.CF01.08.075ISOSAF	8Fr	75cm
ENT.CF01.08.092ISOSAF	8Fr	92cm
ENT.CF01.10.075ISOSAF	10Fr	75cm
ENT.CF01.10.092ISOSAF	10Fr	92cm
ENT.CF01.12.075ISOSAF	12Fr	75cm
ENT.CF01.12.092ISOSAF	12Fr	92cm
ENT.CF01.014.092ISOSAF	14Fr	92cm

PAEDIATRIC FEEDING TUBES




Product Features

- Transparent tube with 2 radiopaque stripes
 - 29-day use
 - Made of polyurethane
- Clearly numbered in centimetres
- Closed rounded tip
 - 2 lateral distal openings
 - ENFit™ male connector with tethered cap
 - 40 units per box

Code	Fr/CH	Length
CA.NCE04PU	4Fr	40cm
CA.NCE05PU	5Fr	50cm
CA.NCE105PU	5Fr	125cm
CA.NCE06PU	6Fr	50cm
CA.NCE106PU	6Fr	125cm
CA.NCE306PU	6Fr	110cm
CA.NCE08PU	8Fr	50cm
CA.NCE108PU	8Fr	125cm
CA.NCE308PU	8Fr	110cm
CA.NCE10PU	10Fr	50cm
CA.NCE110PU	10Fr	125cm
CA.NCE310PU	10Fr	110cm
CA.NCE312PU	12Fr	110cm

LONG TERM NUTRICARE INFANT FEEDING TUBES



Product Features

- Range of Fr sizes (4Fr - 10 Fr)
- Range of lengths (50cm - 125cm)
- Full radiopaque
- Made of polyurethane
- Use up to 90 days
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA free
- 2 non-opposing outlet holes
- 25 units per box

Code	Fr/CH	Length
ENT.NCI04.050ISOSAF	4Fr	50cm
ENT.NCI04.075ISOSAF	4Fr	75cm
ENT.NCI05.050ISOSAF	5Fr	50cm
ENT.NCI05.075ISOSAF	5Fr	75cm
ENT.NCI06.050ISOSAF	6Fr	50cm
ENT.NCI06.075ISOSAF	6Fr	75cm
ENT.NCI06.125ISOSAF	6Fr	125cm
ENT.NCI08.050ISOSAF	8Fr	50cm
ENT.NCI08.075ISOSAF	8Fr	75cm
ENT.NCI08.120ISOSAF	8Fr	120cm
ENT.NCI10.075ISOSAF	10Fr	75cm
ENT.NCI10.120ISOSAF	10Fr	120cm

LONG TERM ADULT FEEDING TUBES



Product Features

- Range of Fr sizes (6Fr - 16Fr)
- Range of lengths (50cm - 120cm) Fully radiopaque
- Made of polyurethane
- Use up to 90 days
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA free
- 2 non-opposing outlet holes
- 5 units per box



Code	Fr/CH	Length	Guide Wire
ENT.NC06.080ISOAF-ROW	6	80cm	Yes
ENT.NC06.110ISOAF-ROW	6	110cm	Yes
ENT.NC08.050ISOAF-ROW	8	50cm	Yes
ENT.NC08.080ISOAF-ROW	8	80cm	Yes
ENT.NC08.092ISOAF-ROW	8	92cm	Yes
ENT.NC08.110ISOAF-ROW	8	110cm	Yes
ENT.NC08.120ISOAF-ROW	8	120cm	Yes
ENT.NC10.092ISOAF-ROW	10	92cm	Yes
ENT.NC10.110ISOAF-ROW	10	110cm	Yes
ENT.NC12.092ISOAF-ROW	12	92cm	Yes
ENT.NC12.110ISOAF-ROW	12	110cm	Yes
ENT.NC12.120ISOAF-ROW	12	120cm	No
ENT.NC14.092ISOAF-ROW	14	92cm	No
ENT.NC14.092ISOAF-ROW-GW	14	92cm	Yes
ENT.NC14.110ISOAF-ROW	14	110cm	Yes
ENT.NC14.120ISOAF-ROW	14	120cm	No
ENT.NC16.092ISOAF-ROW	16	92cm	No
ENT.NC16.092ISOAF-ROW-GW	16	92cm	Yes
ENT.NC16.120ISOAF-ROW	16	120cm	No
ENT.NC08.140ISOAF	8	140cm	Yes
ENT.NC10.140ISOAF	10	140cm	Yes

NUTRICARE FLO FEEDING TUBES - INFANT/PAEDIATRIC POLYURETHANE OPEN END USE



Product Features

- Range of Fr sizes (4 - 12Fr)
- Range of lengths (40cm - 125cm)
- Radiopaque stripe
- Made of polyurethane
- Use up to 30 days
- Clearly numbered in centimetres
- Fixation plaster
- Male ENFit™ connector
- DEHP/Latex/BPA free
- 2 non-opposing outlet holes
- Polyurethane feeding tube with rounded open distal end
- The open end eliminates the collection of feed in the tip of the tube
- Less chance of bacterial colonisation
- 25 units per box



Code	Fr/CH	Length
ENT.NCF04.40	4	40cm
ENT.NCF04.50	4	50cm
ENT.NCF04.75	4	75cm
ENT.NCF05.40	5	40cm
ENT.NCF05.50	5	50cm
ENT.NC05.60	5	60cm
ENT.NC05.75	5	75cm
ENT.NC05.90	5	90cm
ENT.NC05.125	5	125cm
ENT.NC06.40	6	40cm
ENT.NC06.50	6	50cm
ENT.NC06.60	6	60cm

Code	Fr/CH	Length
ENT.NC06.75	6	75cm
ENT.NC06.90	6	90cm
ENT.NC06.125	6	125cm
ENT.NC08.40	8	40cm
ENT.NC08.50	8	50cm
ENT.NC08.60	8	60cm
ENT.NC08.75	8	75cm
ENT.NC08.90	8	90cm
ENT.NC08.125	8	125cm
ENT.NC10.75	10	75cm
ENT.NCF10.110	10	110cm
ENT.NCF12.110	12	110cm

FEEDING TUBES: ADULT LEVIN TUBES - POLYURETHANE



Gastric Drainage & Feeding Tube

Enteral Naso-Gastric Tubes - Levin Type with Guide Wire - Up To 29 Days

Code	Fr/CH	Length
CA.NCE108G	8Fr	110cm
CA.NCE508G	8Fr	130cm
CA.NCE110G	10Fr	110cm
CA.NCE510G	10Fr	130cm
CA.NCE112G	12Fr	110cm
CA.NCE512G	12Fr	130cm
CA.NCE114G	14Fr	110cm
CA.NCE514G	14Fr	130cm
CA.NCE116G	16Fr	110cm
CA.NCE516G	16Fr	130cm

Product Features

- Transparent tube with 2 radiopaque stripes 40% Barium
- Clearly numbered in centimetres
- Closed rounded tip
- ENFit™ male connector with tethered cap
- 4 lateral distal openings, allows for drainage and feeding
- 20 units per box



FEEDING TUBES: ADULT LEVIN TUBES - POLYURETHANE

Gastric Drainage & Feeding Tube




Product Features

- Transparent tube with 2 radiopaque stripes 40% Barium
- Clearly numbered in centimetres
- Closed rounded tip
- ENFit™ male connector with tethered cap
- 4 lateral distal openings, allows for drainage and feeding
- 20 units per box

Enteral Naso-Gastric Tubes - Levin Type - Up To 29 Days

Code	Fr/CH	Length
CA.NCE108AU	8Fr	110cm
CA.NCE508AU	8Fr	130cm
CA.NCE110AU	10Fr	110cm
CA.NCE510AU	10Fr	130cm
CA.NCE112AU	12Fr	110cm
CA.NCE512AU	12Fr	130cm
CA.NCE114AU	14Fr	110cm
CA.NCE514AU	14Fr	130cm
CA.NCE116AU	16Fr	110cm
CA.NCE516AU	16Fr	130cm
CA.NCE118AU	18Fr	110cm
CA.NCE518AU	18Fr	130cm
CA.NCE120AU	20Fr	110cm
CA.NCE520AU	20Fr	130cm

FEEDING TUBES: ADULT DUODENAL/JEJUNAL - POLYURETHANE



Enteral Naso-Duodenal Tubes with Weight and Guide-Wire: Duration of use - Up to 29 days.



Enteral Naso Jejunal Tubes with Weight and Guide-Wire: Duration of use - Up to 29 days



Enteral Naso Duodenal Tubes with Weight: Duration of Use - Up to 29 days

FEEDING TUBES: ADULT DUODENAL/JEJUNAL - POLYURETHANE



Enteral Naso-Duodenal Tubes
with Weight and Guide-Wire:
Duration of use - Up to 29 days.
20 units per box

Code	Fr/CH	Length	Weight
CA.NCE208G	8Fr	110cm	1.1g
CA.NCE210G	10Fr	110cm	1.7g
CA.NCE212G	12Fr	110cm	2.8g
CA.NCE214G	14Fr	110cm	4.0g

Enteral Naso Jejunal Tubes
with Weight and Guide-Wire:
Duration of use - Up to 29 days
20 units per box

Code	Fr/CH	Length	Weight
CA.NCE308G	8Fr	130cm	2.8g
CA.NCE310G	10Fr	130cm	4.0g
CA.NCE312G	12Fr	130cm	4.0g
CA.NCE314G	14Fr	130cm	4.8g

Enteral Naso Duodenal Tubes
with Weight: Duration of Use -
Up to 29 days
20 units per box

Code	Fr/CH	Length	Weight
CA.NCE206AU	6Fr	110cm	1.1g
CA.NCE208AU	8Fr	110cm	1.1g
CA.NCE210AU	10Fr	110cm	1.7g
CA.NCE212AU	12Fr	110cm	2.8g
CA.NCE214AU	14Fr	110cm	4.0g



ENFit™ Sampling Straw (5cm & 15cm)

Stock Code: CA.NCE15S - Wide Bore
Stock Code: CA.NCE05S - Fine Bore

Single Use
150 Per Box



ENFit™ Filter Needle Cap

Stock Code: ENT.FINSO-50

Single Use 5µm 18G x 50mm
30 Per Box



Soft Wide Funnel With Cap

Stock Code: CA.NCE101A

ENFit™ Female
100 Per Box



ENFit™ Stepped Cap

Stock Code: ENT.SBAISO

16-28mm

Stepped Cap for Oral/ENFit™ Syringes
50 per Box



ENFit™ Stepped Cap

Stock Code: CA.NCE100C

10mm to 24mm
100 Per Box



ENFit™ Screw Cap

Stock Code: CA.NCE302C

28mm
50 Per Box



Soft Funnel With Cap

Stock Code: CA.NCE100A

ENFit™ Female

150 per Box



4 - Stepped Connector

Stock Code: CA.CE104A

ENFit™ Male

100 Per Box



5 - Stepped Connector

Stock Code: CA.NCE105A

ENFit™ Male

100 Per Box



Stock code: ENT.BAISO3-S
Single use Size 3: 17-20.5mm
100 per Box

Stock code: ENT.BAISO35-S
Single use Size 4: 19-23mm
100 per Box

Stock code: ENT.BAISO4-S
Single use Size 5: 24-26mm
100 per Box

Internal Bottle Access ENFit™

Allows for tamper proof cap to be replaced

Stock code: ENT.BAISO1-S
Single use Size 1: 12-15mm
100 per Box

Stock code: ENT.BAISO2-S
Single use Size 2: 14-16.5mm
100 Per Box



BOTTLE CAP Enteral Access

Stock Code: CA.NCE400C

Sampling cap with ENFit™
100 Per Box



Y-Port Connector (24 Hour Use)

Stock Code: CA.NCE100Y

Y-Port - 2 Male, 1 Female
100 Per Box



ENFit™ 3-Way Stopcock (24 Hour Use)

Stock Code: CA.NCE300R
3-Way - 2 Male, 1 Female
35 Per Box



ACE Connector (30 Day Use)

Stock Code: DA.485 - Box of 50
Stock Code: DA.485-10 Bag of 10
Ace Connector - Closed Connection System
30 Day Use - Dual Access Port



ENFit™ Drainage Bag (50ml)

Stock Code: ML.DBISO - 50
50ml Drainage Bag
EVA Plastic - Non PVC

20 Per Box



ENFit™ Drainage Bag (500ml) & (1000ml)

Stock Code: MIDMED - 500

500ml Drainage Bag
EVA Plastic - 100cm Tubing - Non PVC
20 Per Box

Stock Code: MIDMED-1000
1000ml Drainage Bag
EVA Plastic - 100cm Tubing - Non PVC
20 Per Box



Gravity Feeding Set with ENFit™ Connector

Stock Code: MIDMED-F81518

PVC Enteral Feeding Bag, DEHP Free, Gravity Use, with Icebag.

1500ml Capacity
Includes ENFit™ side access port
30 Per Box



PH Strips

Stock Code: MN. 92118

Pack of 100



Single Use Cleaning Brush

Stock Code: ENCS2017100

100 Per Box



Multi-Use Cleaning Brush

Stock Code: ENC125212

10 Per Box



Enteral Extension

Stock Code: ENT.ECSF.075-ISOSAF

75cm

25 Per Box

Priming Volume 1.15ml

Stock Code
CA.NCE1315E
2.8mm Ø | 150cm
25 Per Box

Stock Code
CA.NCE1215E
1.6mm Ø | 150mm
50 Per Box

ENTERAL SYRINGE DRIVER



The Enteral Syringe Pump is designed to provide precise and reliable delivery of enteral nutrition and medication. With advanced features and user-friendly design, it ensures optimal patient care and convenience. Featuring a bolus function for controlled, rapid administration. It includes various visual and audible alarms to alert users of important conditions, ensuring safety. Exclusively compatible with Enteral syringes, the pump offers both rate and time delivery modes to suit different medical needs. With a battery life of more than 6 hours at a rate of 5ml/h, it ensures uninterrupted use.

COMPATIBLE WITH:

Our Enteral Syringe Driver can be calibrated to any ENFit™ syringe. Midmed's is currently pre-calibrated to manage the following syringes for your ease

MS.SE05
MS.SE10
MS.SE20
MS.SE60





Features

The MediLime Tummy Button offers a comfortable solution due to its compact size, eliminating the need to carry a tube. We are dedicated to listening to end-users and specialists, continuously innovating our products year after year. The sole purpose of the MediLime Tummy Button is to make living with a feeding system more pleasant and efficient.

- The MediLime Tummy Button consists of pure silicone, is latex- and DEHP free, resulting in healthy skin around the gastro stoma.
- An apple shaped balloon provides an ideal shape for closing the stoma channel.
- Best compatibility with all MediLime Extension Tubes.
- Various balloon diameters (per size 12FR up to 24FR, from 3 up to 8ml volume) and is available in all lengths (cm) (0.8, 1.0, 1.2, 1.5, 1.7, 2.0, 2.3, 2.5, 2.7, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5 and 6.0).
- The Tummy Button comes with an input guide. This can be a useful tool when placing the button if there is scar tissue at the fistula opening or with a larger size of the Tummy Button.



***The Tummy Button is Packaged in a box with these accessories.
Input Guide, Syringe, Lubricating Gel (5 grams), Gauze**

TUMMY BUTTON STOCK CODES



Code	Fr/CH	Size
ML.0320-12-0.8-02	12	0.8cm
ML.0320-12-1.0-02	12	1.0cm
ML.0320-12-1.2-02	12	1.2cm
ML.0320-12-1.7-02	12	1.7cm
ML.0320-12-2.0-02	12	2.0cm
ML.0320-12-2.3-02	12	2.3cm
ML.0320-12-2.5-02	12	2.5cm
ML.0320-12-2.7-02	12	2.7cm
ML.0320-12-3.0-02	12	3.0cm
ML.0320-12-3.5-02	12	3.5cm
ML.0320-12-4.0-02	12	4.0cm



Code	Fr/CH	Size
ML.0320-14-0.8-02	14	0.8cm
ML.0320-14-1.0-02	14	1.0cm
ML.0320-14-1.2-02	14	1.2cm
ML.0320-14-1.7-02	14	1.7cm
ML.0320-14-2.0-02	14	2.0cm
ML.0320-14-2.3-02	14	2.3cm
ML.0320-14-2.5-02	14	2.5cm
ML.0320-14-2.7-02	14	2.7cm
ML.0320-14-3.0-02	14	3.0cm
ML.0320-14-3.5-02	14	3.5cm
ML.0320-14-4.0-02	14	4.0cm
ML.0320-14-4.5-02	14	4.5cm
ML.0320-14-5.0-02	14	5.0cm
ML.0320-14-5.5-02	14	5.5cm
ML.0320-14-6.0-02	14	6.0cm

TUMMY BUTTON STOCK CODES



Code	Fr/CH	Size
ML.0320-16-0.8-02	16	0.8cm
ML.0320-16-1.0-02	16	1.0cm
ML.0320-16-1.2-02	16	1.2cm
ML.0320-16-1.7-02	16	1.7cm
ML.0320-16-2.0-02	16	2.0cm
ML.0320-16-2.3-02	16	2.3cm
ML.0320-16-2.5-02	16	2.5cm
ML.0320-16-2.7-02	16	2.7cm
ML.0320-16-3.0-02	16	3.0cm
ML.0320-16-3.5-02	16	3.5cm
ML.0320-16-4.0-02	16	4.0cm
ML.0320-16-4.5-02	16	4.5cm
ML.0320-16-5.0-02	16	5.0cm
ML.0320-16-5.5-02	16	5.5cm
ML.0320-16-6.0-02	16	6.0cm



Code	Fr/CH	Size
ML.0320-18-0.8-02	18	0.8cm
ML.0320-18-1.0-02	18	1.0cm
ML.0320-18-1.2-02	18	1.2cm
ML.0320-18-1.7-02	18	1.7cm
ML.0320-18-2.0-02	18	2.0cm
ML.0320-18-2.3-02	18	2.3cm
ML.0320-18-2.5-02	18	2.5cm
ML.0320-18-2.7-02	18	2.7cm
ML.0320-18-3.0-02	18	3.0cm
ML.0320-18-3.5-02	18	3.5cm
ML.0320-18-4.0-02	18	4.0cm
ML.0320-18-4.5-02	18	4.5cm
ML.0320-18-5.0-02	18	5.0cm
ML.0320-18-5.5-02	18	5.5cm
ML.0320-18-6.0-02	18	6.0cm

TUMMY BUTTON STOCK CODES



Code	Fr/CH	Size
ML.0320-20-0.8-02	20	0.8cm
ML.0320-20-1.0-02	20	1.0cm
ML.0320-20-1.2-02	20	1.2cm
ML.0320-20-1.7-02	20	1.7cm
ML.0320-20-2.0-02	20	2.0cm
ML.0320-20-2.3-02	20	2.3cm
ML.0320-20-2.5-02	20	2.5cm
ML.0320-20-2.7-02	20	2.7cm
ML.0320-20-3.0-02	20	3.0cm
ML.0320-20-3.5-02	20	3.5cm
ML.0320-20-4.0-02	20	4.0cm
ML.0320-20-4.5-02	20	4.5cm
ML.0320-20-5.0-02	20	5.0cm
ML.0320-20-5.5-02	20	5.5cm
ML.0320-20-6.0-02	20	6.0cm



Code	Fr/CH	Size
ML.0320-24-0.8-02	24	0.8cm
ML.0320-24-1.0-02	24	1.0cm
ML.0320-24-1.2-02	24	1.2cm
ML.0320-24-1.7-02	24	1.7cm
ML.0320-24-2.0-02	24	2.0cm
ML.0320-24-2.3-02	24	2.3cm
ML.0320-24-2.5-02	24	2.5cm
ML.0320-24-2.7-02	24	2.7cm
ML.0320-24-3.0-02	24	3.0cm
ML.0320-24-3.5-02	24	3.5cm
ML.0320-24-4.0-02	24	4.0cm
ML.0320-24-4.5-02	24	4.5cm
ML.0320-24-5.0-02	24	5.0cm
ML.0320-24-5.5-02	24	5.5cm
ML.0320-24-6.0-02	24	6.0cm

STOMA MEASURING DEVICE & GASTRIC MEASURING DEVICE



The MediLime Stoma Measuring Device is used to obtain the correct size of the MediLime Tummy Button. Created from silicone material, the measuring device is PHT and BPA free. As well as latex free, this disposable measuring device is sterile for clinical applications. The stoma canal is measured with this catheter. The use of the Stoma Measuring Device will result in a more precise Tummy Button placement. We strongly recommend using one for each (new) placement of a button in order to optimize the perfect fit.

Stock Code: ML.3025-00



The MediLime Gastric Measuring Device with ENFit™ connection is made for administering enteral nutrition and medication. The Gastric Feeding Tube is made of pure silicone and is therefore skin-friendly and pliable. The size of the tube and the balloon capacity are clearly stated on the shaft.

Code	Fr/CH		
ML.0503-12	12	1-Way	2 Pieces
ML.0503-14	14	1-Way	2 Pieces
ML.0503-16	16	1-Way	2 Pieces
ML.0503-18	18	1-Way	2 Pieces
ML.0503-20	20	1-Way	2 Pieces
ML.0503-22	22	1-Way	2 Pieces
ML.0503-24	24	1-Way	2 Pieces



JEJUNAL FEEDING TUBES & STOMA PLUGS



The MediLime Jejunal Feeding Tube is a 100% silicone device for delivery of enteral nutrition into the distal duodenum or jejunum. The device has a double lumens tube, an air-vented skin disk and a fixation balloon.



Code	Fr/CH			
ML.103080031280A	12	1-Way	2 Pieces	5ml
ML.103080031480A	14	1-Way	2 Pieces	5ml
ML.103080031680A	16	1-Way	2 Pieces	5ml



- Stops stoma sites from closing
- Saves the stoma site until a new gastrostomy tube can be placed
- Includes a range of French gauge sizes in each pack which provides options even when stoma has started to close
- Effective, safe and inexpensive
- Saves staff in emergency situations
- Promptly dealing with gastrostomy emergencies prevents use of other incorrect or inappropriate feeding devices
- Made from high grade medical silicone for flexibility
- Tapered tip for ease of insertion

ENPLUGs are designed for emergency use to keep stomas from closing when a gastrostomy device has fallen out. The ENPLUG packs are available in 2 different lengths and various French sizes in case the stoma has already started to close and the currently used device no longer fits the stoma.

Code	Fr/CH	Size
ENT.EN-PLUGS-S	10,12,14 & 16	4cm
ENT.EN-PLUGS-L	10,12,14 & 16	7cm
ENT.EN-PLUGS-XL	18 & 20	7cm

About Midmed

Welcome to Midmed Pty Ltd, your trusted partner in healthcare excellence! Since 1994, Midmed has been at the forefront of the Australian medical supply industry, delivering top-tier medical equipment and healthcare solutions to businesses and consumers nationwide. With over 25 years of experience, we pride ourselves on our unwavering commitment to quality, innovation, and exceptional customer service.

At Midmed, we understand the unique needs of healthcare facilities, which is why we offer an extensive range of products designed to enhance patient care and streamline operations. From state-of-the-art medical furniture and ergonomic examination tables to reliable patient transfer chairs and advanced medical devices, our comprehensive catalogue is tailored to meet the highest standards of the healthcare industry.

Our dedication to excellence is reflected in every product we offer, ensuring that healthcare professionals can rely on us for the tools they need to provide the best possible care. Join the countless satisfied customers who have made Midmed their go-to source for medical supplies and experience the difference that quality and innovation can make.

Discover the Midmed advantage today and let us help you elevate your healthcare facility to new heights of efficiency and patient satisfaction. Your success is our mission, and we are here to support you every step of the way.





www.midmed.com.au
1300 643 633
sales@midmed.com.au

**Enteral
Feeding**
CATALOGUE 2026