



Montana

Padded, wider armreclinerchair

Designed with wider support arms compared to the Delaware chair, which offers a larger contact surface area. This chair has waterfall back cushions for ultimate comfort, three deeply soft cushions support your head, neck, and full length back.

Suitable for:

- General ward use
- Day care

Fabric in image: Valencia Delft

Specification

	Inches	CM
Seat height to floor	19	48
Seat width	20	51
Seat depth	20	51
Back height	28/29	71/74
Arm height	7	18
Overall height	46	117
Overall width	33	84
Overall width at bottom of arm	28.5	72
Overall depth	34	86
Arm width	6	15

Standard features:

- Waterfall support back design
- Heavy-duty lockable castors
- Available as manual, single motor electric recline only and as a riser recline, with a dual option
- Max user weight 127kg (20stone)

Fabric Options:


- Cadet Shieldplus™
- Valencia
- Haven Shieldplus™
- Dartex (VP)
- Other bespoke fabrics available upon request

Optional Accessories

- Head pillow with weighted tab
- Head roll with weighted tab
- Full length arm cap protective covers
- Accupak rechargeable battery system with alarm (Enable the chair to be wireless - electric chairs only)
- USB charging point to outside arm (will have separate power supply)
- Pressure management seating Liquiform (Med- High) with Dartex VP fabric to seat and leg rest cover
- Wall mounted transformer (avoid having the transformer on the floor)
- Magnet handset (magnetic connection for storage of the handset to outside arm)

Motor Options:

Manual or Electric single Tilt in Space recline only



Single motor Tilt in Space with rise & recline



Dual motor standard with rise & recline

