




**Ultra-low minimum height**  
The industry-leading minimum height of 420 mm enables easier mobilisation & reduced the risk of fall-related injury.



**Electric height adjustment & profiling**  
Fully electric functionality with battery backup enables smooth, comfortable profiling and height adjustment.



**Unrestricted imaging**  
The MultiFlow features a full-length X-Ray tunnel and radiolucent mattress platform for unrestricted imaging.



# MultiFLOW®

## Ultra-Low Height ED Stretcher

Maintain safe patient handling in your emergency department while:

- › Reducing fall-related injuries
- › Enabling easier mobilisation
- › Facilitating faster discharge

Many NHS trusts continue to report delayed discharges and poor patient flow in ED, along with increases in tissue damage and falls due to longer trolley stays.

The MultiFlow stretcher is built to keep patients moving through your hospital, combining an industry-leading low height with a full-length X-Ray tunnel for unrestricted imaging.

### Technical Specification

Safe Working Load	270kg
Max Patient Weight	230kg
Min Platform Height	420 mm
Max Platform Height	820 mm
Platform Dimensions	700 x 2000 mm
ISO 60601-2-52 Compliant?	Yes
X-Ray Compatible?	Yes
Fifth Wheel?	Yes
Configuration Options	X-Ray/Non-X-Ray 2 Section/4 Section platform

Scan for full details

