



PRODUCT COMPARISON

functionalClass

UFSK
INTERNATIONAL

OSYS



implantology
maxillofacial
surgery



ENT surgery



dermatology



500 XLS
Item No.: 1500 0005



I.V.I. treatment



500 ECO
Item No.: 1500 0008

Status:
02/2024



chassis (LxW)

1020 x 710 mm

910 x 600 mm


material

powder-coated welded steel frame

castors

4 castors
2x double castors 100 mm Ø,
360°, braking,
2x 150 mm Ø, fixed,
steering

2x fixed castors, lowerable,
100 mm Ø, rigid,
1x swivel castor 35 mm Ø



braking system

manual

No

central brake

No, brake front castors

No

brake lever

left/right foot lever

No

steering

Yes

No



foot switch

Yes, removable foot control,
standard

Yes, optional:
foot control (Art.Nr.: 3351 1029)



foot switch
functions

on/off (POWER),
adjusting back, seat, height,
3-memory-/ reset-position

2 memory positions
(M1/M2)

top frame

powder-coated
welded steel frame

zinc undercoated, powder-coated
welded steel frame

bed surface (width)

610 mm

600 mm

norm rails

Yes (optional)
Item No.: 3321 0000 / 3320 0000

No

upholstery
(thickness)

back 100 mm
seat 100 mm
leg 80 mm

„skai Tundra“/
„skai Pandoria Plus“

back 110 mm
seat 110 mm
leg 80 mm

headrest

U shape, seams, small hollow

standard: flat

optional: multifunctional headrest
(Item No.: 3351 1034)



headpart adjustments

manual, 1 lever (lengths, tilting)

pushing handles

foot end,
fixed, stainless steel

No





PRODUCT COMPARISON

functionalClass

UFSK
INTERNATIONAL

OSYS



implantology
maxillofacial
surgery



ENT surgery



dermatology



500 XLS

Item No.: 1500 0005



I.V.I. treatment



500 ECO

Item No.: 1500 0008

Status:
02/2024



	Hand control functions	on/off (POWER), adjusting back, seat, foot, height	
	battery display LED	Yes	No
	reset position	Yes	
	memory positions	8	
	CPR button	Yes	
	max. load	200 kg	240 kg
	battery	back up battery Item No.: 3047 0000	No, direct power connection
	mains operation	No	Yes
	surgery cycle	aprox. 100	endless
	entry height	580 mm	550 mm