

HARTMANN



TRAUMA BANDAGE

Your lifesaver
made in Switzerland



Military



Special Forces



Police



Paramedics

The new compression bandage from Switzerland:

Rely on the TRAUMA BANDAGE in emergencies



Fast and effective first aid can save lives. In an emergency, police officers and special forces personnel depend on having the right equipment available to quickly dress wounds and stop bleeding.

In life-threatening situations a quick response can be a matter of life and death. In many cases patients bleed to death because wound care comes too late. The TRAUMA BANDAGE is a genuine lifesaver for special forces personnel and the police. Thanks to its easy handling and effectiveness, the new compression bandage developed by Swiss manufacturer IVF HARTMANN AG can save lives in serious emergencies. It effectively stops heavy bleeding and can even be applied with one hand if needed.

HARTMANN was established in 1818 and has its origins in wound care. The expertise acquired from over 200 years of experience has been incorporated into the development of the TRAUMA BANDAGE. Every part of this compression bandage contributes to the high effectiveness of the product. That applies just as much to the patented pressure applicator as the patented Zetuvit® wound dressing. The TRAUMA BANDAGE ensures fast and reliable wound care and stops heavy bleeding in no time.

Strong hold and optimum fixation

The TRAUMA BANDAGE is available in 3 different widths: 10, 15 and 20 cm. It can be applied to different body parts and is made up of a long-stretch bandage and the trusted Zetuvit® absorbent dressing, which does not stick to the wound thanks to its nonwoven fabric cover. The bandage can be quickly and easily fixed in place with the closure clip. Rows of adhesive dots form an unroll stop to prevent the bandage unrolling. The patented pressure applicator not only effectively stops bleeding but also helps to apply the bandage, meaning that the TRAUMA BANDAGE can be easily and securely applied.

DRESSING

Thanks to the evenly distributed absorbent core and capillary action, blood can be absorbed quickly with the Zetuvit® absorbent dressing. The nonwoven fabric cover prevents the dressing sticking to the wound.



PACKAGING

The vacuum packaging is easily opened with the tear tab and is compact, air-tight, water-tight and gas-tight.



The TRAUMA BANDAGE is part of the basic equipment of the armed forces and emergency services

The TRAUMA BANDAGE is an essential piece of equipment for the armed forces and emergency services. In emergencies with traumatic injuries, it helps prevent patients bleeding out and can save lives.

The predecessors of the TRAUMA BANDAGE have been in use for 100 years. Over decades, the product has been continually refined and improved. The TRAUMA BANDAGE is the result of those many years of optimisation, incorporating the latest findings from research and development.

The new and patented pressure applicator was developed by IVF HARTMANN AG and tested by Empa in St Gallen. Thanks to its flat shape, it distributes pressure evenly, can be perfectly positioned over the wound and helps with looping of the bandage.

PRESSURE APPLICATOR

The newly developed pressure applicator ensures even distribution of pressure and optimum positioning over the wound. It also simplifies looping of the bandage.

CLOSURE CLIP

The bandage can easily be fixed in place with the closure clip.

UNROLL STOP

The unroll stop prevents the bandage from inadvertently unrolling. This ensures hygienic application.

BANDAGE

The long-stretch bandage contains 65% cotton and can deliberately generate pressure.

The TRAUMA BANDAGE

- ✓ The components are a bandage, the Zetuvit® absorbent dressing, a pressure applicator, an unroll stop and a closure clip.
- ✓ The large tear tab makes it easy to open the compact packaging, even when wearing gloves.
- ✓ Pressure can be generated with control using the long-stretch bandage.
- ✓ The pressure applicator ensures not only even distribution of pressure but also easy and efficient application.
- ✓ The Zetuvit® absorbent dressing does not stick to the wound.
- ✓ The TRAUMA BANDAGE is a sterile class Is (MDR) medical device that is simply packed for quick application. The packaging is airtight, water-tight and gas-tight.

Place your order now!

Your contact for enquiries:

DOMMENIK STEPPUTAT

+41 79 928 01 40

dommenik.stepputat@hartmann.info

Your premium product in Swiss Made quality

Always operational with the TRAUMA BANDAGE

The TRAUMA BANDAGE is the perfect choice for all those who need a haemostasis solution that is both highly effective and available at all times. The TRAUMA BANDAGE covers wounds, protects them against soiling and effectively stops bleeding.

The TRAUMA BANDAGE is available not only as an individual product but also combined with a tourniquet and single-use gloves in the practical Trauma-Box® bag or Trauma-Box® dispenser.

Anyone can save a life with the Trauma-Box®

Worldwide, more people die from the consequences of wounds or injuries not being treated or being treated too late than from diseases such as malaria, tuberculosis or HIV¹. One of the most common, and preventable, causes of death is blood loss.

If large, open wounds are not attended to promptly, the patient can bleed out in just a few seconds. That is why stopping critical bleeding can be a matter of life and death. The Trauma-Box® means first-aiders can intervene directly and prevent serious consequences in emergencies. IVF HARTMANN AG has developed the Trauma-Box® in collaboration with experts from the German Trauma Foundation. The products in the Trauma-Box® help to stop bleeding until medical professionals arrive.

¹ German Trauma Foundation brochure

All in one: the Trauma-Box®

Emergency services personnel can keep the Trauma-Box® bag in their service vehicle, carry it with them or have it permanently available in the Trauma-Box® dispenser in high-risk areas. The Trauma-Box® dispenser can be placed prominently in locations where large numbers of people gather, such as reception areas, railway stations, shooting ranges or shopping centres. As a first aid kit for emergency haemostasis, the practical Trauma-Box® bag in the dispenser is ready for use at all times.

In an emergency, employees and passers-by can quickly access the Trauma-Box® and the TRAUMA BANDAGE and provide first aid at any time. The eye-catching Trauma-Box® dispenser is weatherproof and vandal-proof and is quick and easy to open in an emergency.

Contents of the Trauma-Box®

- ✓ TRAUMA BANDAGE
- ✓ Tourniquet
- ✓ Single-use gloves



HARTMANN



Helps. Cares. Protects.

IVF HARTMANN AG
Victor-von-Bruns-Strasse 28
Postfach 634
CH-8212 Neuhausen

Tel. +41 52 674 31 11
Fax +41 52 672 74 41
info@ivf.hartmann.info
www.ivf.hartmann.info