

CO.E200IX TRAUMA STRETCHER



www.midmed.com.au




1300 643 633

CO.E200IX TRAUMA STRETCHER

Cogent's Trauma Stretcher with added knee gatch position & 4 section backboard, is very suitable for use in different transport situations. The standard platform provides trauma treatment to severe patients, and can efficiently transport patients. The full-length X-ray allows the X-ray cassette tray to slide from head to foot sections to do X-ray examinations, offering flexibly without moving patients.





STANDARD FEATURES

- Integrated storage tray with Oxygen tank through which is made with ABS plastic
- Foot pedal for hydraulic system control on both sides
- Structure with four-section bed board
- Drainage Bag clasp x 2 & Urinary Bag clasp x2
- Foot controlled brake / steer pedal for four wheels at two ends of stretcher
- IV pole inserted hole x 5
- Dual pneumatic assisted backrest
- Two stage movable IV pole x 1
- Two stage permanent folding IV pole x 1
- Foot-end assist push bar
- Velcro tape underneath the mattress
- Protective corner bumpers x 4
- Hook / Jack for upright oxygen tank holder
- Retractable steering fifth wheel
- Full-length X-ray cassette slide tray
- No-Crank pneumatic knee gatch
- Upright Oxygen tank holder
- Collapsible siderail (PP plastic)
- Cylinder sleeves (PP plastic)
- Head-end ergonomic push handles

SPECIFICATIONS

- Overall Length: 213cm
- Overall Width: 85cm
- Height Range (from ground to litter top):
High: 99cm ± 2
Low: 64cm ± 2
- Weight Capacity: 318kg
- Litter Positioning:
Backrest: 0-90° ± 2°
Trend. / Reverse Trend: 0-18° ± 2°
- Mattress Base: 195 x 64cm
- Mattress Dimensions:
195 x 72 x 10cm (blue)
- Siderail: 126 x 39cm
- Caster Diameter: 20cm
- Manual Knee Gatch 0-35° ± 2°
- Two years warranty for stretcher under normal use after date of shipping

Ordering Code: CO.E200IX



www.midmed.com.au



1300 643 633