

Marina automatic recliner



The *Marina* Automatic recliner is an advancement on DECAM's original product which has demonstrated outstanding reliability and durability through use by hundreds of customers in both private and public hospital settings.

The Marina Automatic recliner embodies a novel concept, namely an automatic reclining chair for indoor use which is also particularly useful for patient mobility.

○ ● ● ● ● positions

reclined position

thanks to two independent motors, the Marina Automatic recliner features multiple reclining positions controlled with an easy to use, hand held controller



upright position

the armrests can be raised up to transfer the patient to or from a bed without any effort. The armrests are movable in any chair position



sleeper

with a simple touch of a button and without any effort, the Marina Automatic recliner can be converted into a bed which provides an exceptional level of comfort



Lifting position (optional)

Patients often suffer from reduced mobility. Helping patients get out of the chair will not only make life easier for the patients but will also reduce caregiver injuries



○ ● ● ● ● standard features

push bar

an innovative design that facilitates the movement and rotation of the chair

urethane arm caps

provides users with a comfortable grip and protection for their arms

moisture barrier

offers protection from fluids such as blood, sweat and urine, with minimum maintenance..

multiple accessories bar

allows for the attaching many different accessories such as an IV pole, a blood donor arm support, hanging or drainage bags or an articulating table

controller

easy to use, hand held controller

moveable armrests

the armrests can be lifted to transfer the patient to or from a bed without any effort



push bar and arm cap palette
personalize your chair with wengue, orange or anthracite colours



pull out food rest
for additional support for legs when seated

zero turning radius
eases mobility in small spaces

casters
silent and resistant double casters with individual locks

○ ● ● ● ● optional features

battery

a rechargeable battery provides the user with a much greater level of autonomy

oxygen tank holder

attaches to the back of the chair to carry an oxygen cylinder

swing away table

allows users with limited mobility with greater flexibility to eat, read or compute

head rest

the adjustable head rest is removable, and provides greater comfort to the user

blood donor arm holder

the adjustable donor arm holder provides a comfortable arm position for a donor or patient

IV pole

stainless steel pole that fits into the adjustable IV pole holder

IV Pole holder

attaches to the multiple accessories bar, and is height, tilt and direction adjustable

○ ● ● ● ● technical specifications

Marina automatic

dimensions	overall		seat		arm height		full recline	
width	62,5 cm	24,4"	50 cm	19,7"	72 cm	28,3"		
depth	80 cm	31,5"	49 cm	19,3"			170 cm	66,9"
height	128 cm	50,4"	52 cm	20,5"				
weight	69Kg	151,8 lb						
max load	150 Kg	330 lb						

Marina automatic (wider seat)

dimensions	overall		seat		arm height		full recline	
width	67 cm	26,4"	56 cm	22"	72 cm	28,3"		
depth	75 cm	29,5"	49 cm	19,3"			170 cm	66,9"
height	127 cm	50"	52 cm	20,5"				
weight	75 Kg	165 lb						
max load	150 Kg	330 lb						

○ ● ● ● ● fabrics

The Decam textile collection, based on Spradling® Diamante and Silvertex patterns is designed specifically to meet the strictest requirements of the healthcare environment by creating **effectives barriers against germs, abrasion and stains**



* visit our website for a full list of available colors

Our standard collection meets the following features:

- | | | | |
|--------------------------|------------------|------------------|--------------------------|
| Abrasion resistance | Urine protection | Sweat Protection | 100% Polyester |
| Antibacterial Protection | Blood Protection | Mildew Resistant | Sulphide Stain Resistant |
| Antifungal Protection | Stain Protection | | |

This brochure provides an overview of the main and benefits of the DECAM product. For a full list of products available, including specifications, fabrics, and finishes, please visit www.decamlux.com. For pricing enquiries contact your local DECAM sales representative.